



Soleo*line*

Modern technologies
for electro-, ultrasound-
and combination therapy

Zimmer
MedizinSysteme

NEW GENERATION – PROVEN THERAPY PRINCIPLES

THE NEW SOLEO SERIES

What is the Soleoline?

A state-of-the-art, innovative product line with three different device variants.

SOLEO GALVA

A state-of-the-art, innovative electrotherapy device with the possibility of additional use of the Vaco-S.

SOLEO SONO

A state-of-the-art, innovative ultrasound therapy device.

SOLEO SONOSTIM

A state-of-the-art, innovative combination device for electrotherapy and ultrasound therapy with the possibility of additional use of the Vaco-S.





4 DEVICES – 4 REASONS TO BUY



COMPLETE THERAPY UNIT

- » Adapted to international standard specifications
- » Fully independent therapy channels usable for pain management,
- » Suitable for pain management, muscle stimulation, and muscle strengthening
- » Ultrasound therapy 1 & 3 MHz
- » Combination therapy between ultrasound and electrostimulation



FASTER

- » New microcontroller, new touch
- » Clear, clear menu structure (graphical user interface)
- » All relevant and important information on one screen
- » Integrated protocols and various help screens: Avoid losing valuable treatment time!
- » The therapist can concentrate on the patient, while the 4 series takes care of the treatment!



THE RIGHT DECISION WITH RESPONSIBILITY

- » Integrated and evidence-based protocols
- » Adapt to the needs of customers and patients
- » A step toward effective and reliable treatment!
- » Helps you make a rational therapeutic decision!



EASY TO USE

- » Large LCD screen with touch
- » Modernized and simplified user interface
- » Battery operation as an option

PRECISION & VERSATILITY – IN ANY FORM OF THERAPY

••• Soleo *Galva*

Versatile electrotherapy for rehabilitation and sports

Wide range of applications

Soleo Galva offers a wide range of stimulation current and special current forms – ideal for individual therapy approaches in medical rehabilitation. Soleo Galva – modern electrotherapy for the highest demands.

FUNCTIONS

» EXTENSIVE SELECTION OF CURRENT WAVEFORMS:

Interference current, premodulated current,
TENS current, dia-dynamic current, galvanic current,
Mirco current, high voltage, IG current

» 2 COMPLETELY INDEPENDENT THERAPY CHANNELS

Twinstim: 1 program on 2 channels
Duostim: 2 programs simultaneously
For flexible and time-saving treatments.

» LARGE 10 INCH LCD SCREEN WITH TOUCH

» MODERNIZED AND SIMPLIFIED USER INTERFACE

» OPTIONAL BATTERY USABLE WITHOUT CONNECTION

» CAN BE COMBINED WITH THE VACO-S



Soleo *Sono*

Ultrasound therapy rethought

Effective in a targeted manner – where the pain arises

Soleo Sono offers modern ultrasound technology for effective and precise therapy – ergonomic, reliable and versatile for acute and chronic complaints.

Fast, secure and intuitive

- » Intuitive therapy selection by indication
- » More than 50 treatment options with recommendations
- » 120 memory locations for favorites and individual programs

Efficient in application

With its sophisticated operation and wide range of indications, Soleo Sono is a powerful ultrasound device for demanding everyday practice.

FUNCTIONS

- » MULTI-FREQUENCY ULTRASOUND 1 & 3
- » MHZ 2 CONNECTIONS FOR THE HANDPIECES
- » LARGE 10 INCH LCD SCREEN WITH TOUCH
- » MODERNIZED AND SIMPLIFIED USER INTERFACE
- » OPTIONAL BATTERY USABLE WITHOUT CONNECTION



ULTRASOUND & ELECTRICITY – PERFECTLY COMBINED

••• Soleo *SonoStim*

Ultrasound and stimulation current in perfect combination

Maximum versatility in one system

Soleo SonoStim combines the comprehensive functions of Soleo Sono (ultrasound therapy) and Soleo Galva (electrotherapy) in a compact, mobile, high-tech device – including the combination therapy of ultrasound and stimulation current. With over 300 protocols for various indications, it offers an impressive range of therapeutic options.

Simultaneous therapy – double effect, one goal

Thanks to the simultaneous delivery of ultrasound and stimulation current, the therapeutic effects of both forms of energy can be specifically combined. Simultaneous therapy is particularly effective for therapy-resistant complaints and functional disorders – for even faster and more sustainable therapeutic success.

FUNKTIONEN

» 3 COMPLETELY INDEPENDENT THERAPY CHANNELS

» EXTENSIVE SELECTION OF CURRENT WAVEFORMS

Interference current, premodulated current, TENS current, dia-dynamic current, galvanic current, Mirco current, high voltage, IG current

» MULTI-FREQUENCY ULTRASOUND 1 & 3 MHZ

» 2 CONNECTIONS FOR THE HANDPIECES

» LARGE 10 INCH LCD SCREEN WITH TOUCH

» COMBINATION THERAPY

» OPTIONAL BATTERY USABLE WITHOUT CONNECTION

» CAN BE COMBINED WITH THE VACO-S



Soleo *Vaco-S*

Therapy with comfort

The combination of electrotherapy and a pleasant suction massage allows a comfortable and effective application.

FUNKTIONEN

- » 2-CHANNEL VACUUM UNIT
- » CONTINUOUS OR PULSED MODE
- » BUILT-IN WATER RESERVOIR
- » CONNECTORS ON THE REAR



FROM ACUTE TO CHRONIC – THERAPY AS NEEDED

••• TREATMENT GOALS

» PAIN MANAGEMENT

Acute
Subacute
Chronic

» MUSCLE STIMULATION

» MUSCLE STRENGTHENING

» IMPROVEMENT OF BLOOD CIRCULATION

MUSCLE STIMULATION

MUSCLE
STRENGTHENING

IMPROVEMENT OF
BLOOD CIRCULATION

PAIN RELIEF

••• SAFETY AND USAGE INSTRUCTIONS::

- » Regular maintenance and hygienic cleaning required.
- » Read the operating instructions before use.
- » Only use on intact skin.
- » This is a medical device. It complies with all applicable medical regulations.

••• WARNINGS & RISKS::

- » Caution in cases of hemorrhage, circulatory disorders, loss of sensation
- » Risk of blistering, aseptic necrosis, eye damage
- » Unclear effects on thyroid gland & lymph nodes
- » Only use applicator on intact skin
- » Do not apply directly to wounds, only to the edges
- » Ultrasound should not be used on unconscious patients



SCOPE OF DELIVERY

1X SOLEO GALVA WITH:

1x power cable, 2x electrode cable, 2x alligator clips red, 2x alligator clips black, 2x rubber electrodes 2x sponge pocket, 1x mixed single-use electrodes, 2x velcro tape

1X SOLEO SONO WITH

1x Ultrasound head 1 and 3 MHz, 5 cm², 1x power cable, 1x holder for ultrasound head, 1x bottle Sono gel 0,5l

1X SOLEO SONOSTIM WITH:

1x ultrasound head 1 and 3 MHz, 5 cm², 1x power cable, 2x electrode cable, 2x alligator clips red, 2x alligator clips black, 2x rubber electrodes, 2x sponge pocket, 1x mixed single-use electrodes, 2x velcro tape, 1x holder for ultrasound head, 1x bottle Sono gel 0,5l

OPTIONAL VACO-S DEVICE:

1x connecting cable, 2x electrode sleeves black, 2x electrode sleeves red, 2x packs Vaco electrodes, 2x packs Viscose sponge



OPTIONAL

ACCESSORIES

Adhesive electrodes small and large, Sono gel , battery

TECHNICAL DETAILS –

COMPACTLY

TECHNICAL DATA

Ultrasound therapy

Ultrasound

Applicators

Frequencies

Small applicator

Maximum power output

Intensity levels

Large applicator

Maximum power output

Intensity levels

Protocol

Sono

1-channel ultrasound

1 cm² und 5 cm²

1 MHz und 3 MHz

1 cm², ERA = 1 cm²

2 W / cm²

0–3 W / cm²

5 cm², ERA = 4 cm²

8 W / cm²

0–3 W / cm²

Over 50 treatment recommendations



Electrotherapy

Number of channels

Protocol

Galva

Two, completely independent (intensity and parameters)

Over 120 treatment recommendations

Combination therapy

Protocol

SonoStim = Sono + Galva

Over 300 treatment recommendations

Soleoline in general

Favorite memory

120 entries in the favorites list

Dimensions / weight

325 mm (L) x 242 mm (W) x 192 mm (H) $\pm$ 5 mm / 2.9 kg

Vacuum unit

Possible applications

continuous / pulsed

Control

Via the touch screen

Dimensions / weight

356 mm (L) x 227.5 mm (W) x 84 mm (H) $\pm$ 5 mm / 1.8 kg



Zimmer

MedizinSysteme

Zimmer MedizinSysteme GmbH
Junkersstraße 9
89231 Neu-Ulm, Germany
Tel. +49 7 31. 97 61-291
Fax +49 7 31. 97 61-299
export@zimmer.de
www.zimmer.de



MD

EN | NM | 10105699 | 0825 | Version 1 | Subject to change

